

Differential Display

Made Simple with

QuickPoint™ Pre-Cast Gels!

NOVEX makes differential display analysis a breeze by introducing the QuickPoint™ Rapid Nucleic Acid Separation System:

- **SMALL**
10 x 12cm gels
- **FAST**
20 minute runs
- **PRE-CAST!**
6 month shelf-life

The QuickPoint™ system provides a fast, easy-to-use, pre-cast gel capable of resolving differentially displayed products in less than 20 minutes! Just dilute your display products in QuickPoint Sample Loading Buffer, load onto a QuickPoint pre-cast gel, run, and quantitate by autoradiography. QuickPoint gels provide fast and accurate differential display results every time without the need to pour your own cumbersome sequencing gel!

GenHunter's RNAimage™ Kit utilizes primers designed to reverse transcribe the 5' poly A tail portions of mRNA's, followed by PCR amplification with additional upstream arbitrary primers. The mRNA sub-populations are then visualized by denaturing polyacrylamide gel electrophoresis. This enables direct side-by-side comparison of the mRNAs between two or more related cells and allows differentially expressed genes to be identified.

By utilizing QuickPoint Pre-Cast Gels along with the RNAimage Kit, differential display is fast, easy and convenient.

Comparison of 2 oligo-dT primers from the GenHunter RNAimage™ Kit.

The RNA Differential Display process is covered by U.S. Patent 5,262,311 and pending foreign patents, licensed to GenHunter Corporation.
PCR process is covered by patents owned by Hoffman-LaRoche, Inc. RNAimage™ is a registered trademark of GenHunter Corporation. For information relating to the RNAimage™ Kit, please contact GenHunter at 1-800-311-0260.



**Differential
Display**

Made Simple with

QuickPoint™ Pre-Cast Gels!



United States

NOVEX

11040 Roselle Street

San Diego, CA 92121

Tel: (619) 452-6634 / (800) 456-6839

Fax: (619) 452-6635

Technical Service: (800) 556-6839

E-mail Address: nvxinfo@novex.com

Germany

NOVEX Electrophoresis GmbH

Brüningstrasse 50

Gebäude C584

D-65929 Frankfurt/Main, Germany

Tel: +49 (69) 3058 2100 / +800 2626 3333

Fax: +49 (69) 3058 2102

E-mail Address: info@novex.de

UK Technical Service: 0800 965822

UK Sales: 0800 965823

UK Fax: 0800 965824

Website: <http://www.novex.com>

NOVEX, NOVEX, **QuickPoint**, QuickPoint, are trademarks of Novel Experimental Technology. ©1998 Novel Experimental Technology. All Rights Reserved.

1-800-456-6839

EP066 10/98

PROTOCOL

PLEASE READ THE QUICKPOINT GEL INSTRUCTION BOOKLET BEFORE PROCEEDING.

1. Pre-run the QuickPoint gel for 5 minutes according to the recommended protocol.
2. Mix the PCR reaction (20µl) with an equal volume (20µl) of the QuickPoint Sample Loading Buffer.
3. Heat samples at 80°C for 2 minutes immediately before loading onto the pre-run QuickPoint gel.
4. Load 1µl of the prepared samples onto the QuickPoint gel. To assure a straight banding pattern, 0.5µl of the QuickPoint Sample Loading Buffer should be loaded in remaining empty wells.
5. Run gel until the bromophenol blue dye reaches the bottom of the gel. This should take approximately 10 minutes. Longer run times may be necessary depending upon fragment size.
6. Remove the glass cassette from the cell. Using a razor blade, slice tape and open cassette. The gel will remain on the notched glass plate.
7. Wash the gel/glass plate in water for 5 minutes on a low speed shaker. Do NOT fix the gel.
8. Dry the gel in an 80°C oven for 20 minutes.
9. Expose, gel side up to X-ray film overnight.
10. Bands of interest can be reamplified. NOVEX recommends using the Micropure 0.22 separators and Gel Nebulizers™ from Millipore (Cat# 42600) with the following protocol to extract the DNA band from the QuickPoint gel.
 - a. Orient the autoradiogram with the gel using the dark outline of the wells as a guideline.
 - b. Rehydrate the gel by covering just the band of interest with 10µl H₂O.
 - c. Cut and scrape out the rehydrated band of interest with a spatula and place into the bottom of a Gel Nebulizer™; then into the Micropure 0.22 separator.
 - d. Add 10µl of H₂O on top of the gel slice. Soak the gel slice for 5 minutes.
 - e. Microcentrifuge at 14,000rpm for 10 minutes.
 - f. Approximately 10µl of extracted DNA sample will collect in the vial. Use 4µl for reamplification.

Ordering Information

Product

QuickPoint™ Starter Kit

Kit includes: 1 box of 4 QuickPoint Gels (QP9731),

QP Running Buffer, 50X, 250ml (QP9732),

QP Sample Loading Buffer, 2.5ml (QP9733).

QuickPoint™ Cell

QuickPoint Gels (16 well, 250µm thick, 4/box)

QuickPoint Running Buffer, 50X, 250ml

QuickPoint Sample Loading Buffer, 2.5ml

QuickPoint Kit

Kit includes: 2 boxes of 4 QuickPoint Gels (QP9731),

QP Running Buffer, 50X, 250ml (QP9732), QP Sample Loading Buffer, 2.5ml (QP9733).

Code

QP9735

E19700

QP9731

QP9732

QP9733

QP9736