


E-Gel™ 96—Automated agarose gel electrophoresis as easy as Plug & Play




The E-Gel® 96 system is:

- **Ultra-fast** — Analyze 96 samples in just 12 minutes
- **Efficient** — Run multiple gels simultaneously, over 20,000 samples per day
- **Convenient** — Compatible with 8-, 12-, and 96-tip robot loading for automated operation



To experience the latest in high-throughput agarose gel electrophoresis, check out the E-Gel® 96 system on the enclosed CD.



Invitrogen
life technologies
1600 Faraday Avenue
Carlsbad, CA 92008

PRSR1 STD
U.S. Postage:
PAID
Permit No. 825
San Diego, CA

Plug & Play electrophoresis

For high-throughput electrophoresis, simply place the pre-cast E-Gel[®] 96 cassette on an E-Gel[®] 96 base (A), load your samples (B), and run (C). In just 12 minutes, you'll get results (D).

A. Set up



B. Load samples



C. Run



D. Get results



To learn how you can turn your time-consuming, labor-intensive agarose gel electrophoresis assignments into an automated, high-throughput process, check out the E-Gel[®] 96 system on the enclosed CD.

For more information on E-Gels[®] and the E-Gel[®] 96 system, visit our web site at www.invitrogen.com/egels.

Get plugged in today

Streamline your high-throughput DNA agarose electrophoresis applications. Order the E-Gel[®] 96 system today.

Product	Quantity	Cat. no.
E-Gel [®] 96 1% agarose gel	8 gels	G7008-01
E-Gel [®] 96 2% agarose gel	8 gels	G7008-02
E-Gel [®] 96 mother base	1 base	G7100-01
E-Gel [®] 96 daughter base	1 base	G7200-01
E-Gel [®] 96 holder	2 holders	G7300-01
E-Gel [®] 96 Editor software	1 copy	*FREE



*The E-Gel[®] 96 Editor software is available free of charge with purchase of E-Gel[®] 96 gels and related equipment. This software may be downloaded from the Invitrogen web site at www.invitrogen.com/egels.

E-Gel[®] 96 products are covered by Limited Label License No. 61. Please refer to the Invitrogen web site or catalog for the corresponding Limited Label License.

Toll Free: 800 955 6288